


HEGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°5/6



Ref. IS01405



30/07/2025 06:20:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

HEGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°5/6



Ref. IS01405

Description

Hegar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029566	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.2
Packed unit / Weight	0.043	Bare Unit / Weight (kg)	0.043
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/07/2025 06:20:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr