

HEGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°13/14



Ref. IS01413



29/07/2025 23:11:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

HEGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°13/14



Ref. IS01413

Description

Hegar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029603	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.3
Packed unit / Weight	0.063	Bare Unit / Weight (kg)	0.063
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/07/2025 23:11:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr