

KOENING LIGATURE CONDUCTOR 19.5 CM X 3 MM

Ref. IS02615



30/07/2025 09:31:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

KOENING LIGATURE CONDUCTOR 19.5 CM X 3 MM



Ref. IS02615

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040087511	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	6.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23.6
Packed unit / Weight	0.042	Bare Unit / Weight (kg)	0.042
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 09:31:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr