

PINARD OBSTETRICAL STETHOSCOPE, ALUMINIUM

Ref. IS07601



30/07/2025 03:59:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PINARD OBSTETRICAL STETHOSCOPE, ALUMINIUM



Ref. IS07601

Stethoscope foetal Pinard used to listen to foetal heart sounds

Description

Small conical instrument, of which the top allows to the practitioner to place the ear, while the flat and circular basis is placed on the abdomen of the pregnant woman.

Designed for frequent and disinfection with hospital-grade products.

- Single chestpiece
- Lightweight
- Unbreakable aluminum
- Easy to handle
- Easy to clean

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030548	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	7
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	14
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 03:59:25 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr