

## DOYEN NEEDLE HOLDER 15 CM

Ref. IP30415



30/07/2025 08:38:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DOYEN NEEDLE HOLDER 15 CM



Ref. IP30415

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 3661040024837 | HS Code                   | 90189084 |
| Warranty (months)           | 24            | Country                   | Pakistan |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Bare Unit / Length (cm)   | 15       |
| Packed unit / Weight        | 0.039         | Bare Unit / Weight (kg)   | 0.039    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Packed unit / Quantity      | 1             |                           |          |

30/07/2025 08:38:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)