

BALLENGER SPONGE FORCEPS STRAIGHT 18 CM (STANDARD EYE)



Ref. IP45618



30/07/2025 09:25:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

BALLENGER SPONGE FORCEPS STRAIGHT 18 CM (STANDARD EYE)



Ref. IP45618

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040025858	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.7
Packed unit / Weight	0.046	Bare Unit / Weight (kg)	0.046
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/07/2025 09:25:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr