

DUJARDIN EYE BALL FIXATION FORCEPS, 11.5CM

Ref. IP52100



29/07/2025 23:46:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DUJARDIN EYE BALL FIXATION FORCEPS, 11.5CM



Ref. IP52100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040027784	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	7.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20
Packed unit / Weight	0.013	Bare Unit / Weight (kg)	0.013
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:46:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr