

DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 5 CM STRAIGHT

Ref. IP54005



29/07/2025 23:29:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 5 CM STRAIGHT



Ref. IP54005

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040028194	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	5.5
Packed unit / Weight	0.008	Bare Unit / Weight (kg)	0.008
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/07/2025 23:29:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr