

SEQUESTRUM FORCEPS STRAIGHT 23 CM

Ref. IP54423



 30/07/2025 01:15:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEQUESTRUM FORCEPS STRAIGHT 23 CM



Ref. IP54423

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040028316	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23
Packed unit / Weight	0.226	Bare Unit / Weight (kg)	0.223
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 01:15:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr