

BAR UMBILICAL CLAMP CUTTING FORCEPS, SEMI-SHARPED

Ref. IP56401



30/07/2025 05:40:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

BAR UMBILICAL CLAMP CUTTING FORCEPS, SEMI-SHARPED



Ref. IP56401

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040087016	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	3.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19.2
Packed unit / Weight	0.112	Bare Unit / Weight (kg)	0.112
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/07/2025 05:40:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr