

DOYEN ADULT RASPIRATORY RIGHT SIDE 17 CM

Ref. IR06400



29/07/2025 23:55:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

DOYEN ADULT RASPIATORY RIGHT SIDE 17 CM



Ref. IR06400

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029078	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	3	Bare Unit / Width (cm)	3
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.4
Packed unit / Length (cm)	17	Bare Unit / Length (cm)	17
Packed unit / Weight	0.07	Bare Unit / Weight (kg)	0.07
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/07/2025 23:55:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr