

## THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

Ref. DWPRO6K08

**WelchAllyn®**

Advancing Frontline Care™

CE  
0297

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 YEARS  
GUARANTEE



01/08/2025 22:42:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

**WelchAllyn®**

Advancing Frontline Care™

Ref. DWPRO6K08

**Tympanic infrared thermometer for pros, in 2 seconds and with Electric Charging Station**

## Description

Welch Allyn's Braun ThermoScan PRO 6000 ear thermometer, featuring the most advanced measurement technology, provides the speed and simplicity healthcare professionals need to capture accurate, reliable and comfortable temperature readings for patients.

- Measurement time: 2 - 3 seconds

- The innovative PerfectTemp™ technology improves the accuracy of measurements by detecting the orientation and depth of insertion of the probe into the ear canal and adjusts the temperature calculation parameters.

- ExacTemp™ technology contributes to the accuracy of the measurement by detecting the stability of the probe during the measurement

- It has a 60-second pulse counter for manual measurement of pulse rate and breathing, stopwatch, buttons to recall the last measurement and change of graduation C/F

- Type of display: Large backlit liquid crystal display, four digits plus specific icons

- Automatic standby: About 10 seconds after the last measurement

- Power: rechargeable battery thanks to its charging station).

- Probe covers for ThermoScan® Pro 6000 protect the probe and help minimize the risk of cross-contamination.

- Small: 152 mm × 44 mm × 33 mm (6" × 1.7" × 1.3") and lightweight: 100 grams (without batteries).

- Battery life: 1000 measurements (AA alkaline batteries)

- Comes with: short base, charging station, rechargeable battery and a box of 20 probe covers

 01/08/2025 22:42:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## THERMOSCAN PRO 6000 TYMPANIC THERMOMETER WITH STATION

**WelchAllyn®**

Advancing Frontline Care™

Ref. DWPRO6K08

### Features

CE Marking	CE0297	Medical device	Class IIa
Bar code	3661041009826	HS Code	90251900
Warranty (months)	36		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Bare Unit / Width (cm)	3.3
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	15.2
Packed unit / Length (cm)	25	Bare Unit / Length (cm)	4.4
Packed unit / Weight	0.46	Bare Unit / Weight (kg)	0.44
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

01/08/2025 22:42:43 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)