

GP DOCTOR BAG BLACK

Ref. MEGP03_R

 **ELITE BAGS**



 01/08/2025 23:01:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

GP DOCTOR BAG BLACK



Ref. MEGP03_R

Description

Lightweight and easy to use bag thanks to its internal velcro dividers that allow to organize the sanitary equipment in an easy and fast way. Its system of integrated shoulder straps makes it possible to transform it into a backpack.

-Light and resistant textile material, anti-tear (Ripstop)

- Special area for personal belongings, business cards, pen, etc ...

-Document compartment, A4

-Internal velcro dividers

- Insulated amp maker, 43 bulbs

- Zipper with eyelet for padlock.

- Large external pocket with double compartment and elastic bands

-Integrated shoulder straps for transformation into a backpack

-Weight: 1 kg

-Outdoor measures: 40 x 21 x 25 cm / 21 L

-Maximum recommended capacity: 3.5 Kg

- Sold empty

- Available in black (MEGP02_N), anthracite/salmon (MEGP03_R) and grey/turquoise (MEGP03_B).

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540461	HS Code	42029900
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	21
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	25
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	40
Packed unit / Weight	1.5	Bare Unit / Weight (kg)	1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Anthracite/saumon		

01/08/2025 23:01:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr