

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330012

05/08/2025 22:55:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330012

Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

- Adult (M) cotton, velcro, grey: 330011
- Adult (M) cotton, velcro, red : 330013
- Adult (M) navy, with straps: 330112
- Adult (M) red, with straps: 330113
- Adult (L) cotton, velcro, navy: 330002
- Child/small adult (S), cotton, velcro, grey: 330021

Features

CE Marking	NONE	Medical device	No
Bar code	3700446028541	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.19
Bare Unit / Weight (kg)	0.19		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


05/08/2025 22:55:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr


+33(0)4 42 90 31 31

info@gsh-med.fr