

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513229



CE
0459

DISPOSITIF MEDICAL
CLASSE
I_m
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

12/08/2025 22:25:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)



Ref. 513229

Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007195	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	7
Packed unit / Length (cm)	19	Packed unit / Weight	0.526
Bare Unit / Weight (kg)	0.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

 12/08/2025 22:25:11 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr