

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)

Ref. 520627



CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**I<sub>m</sub>**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 YEARS  
GUARANTEE

12/08/2025 21:26:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LIAN NM / MULTI CUFFS / NYLON / (L, M, S)



Ref. 520627

=====

### Description

#### HIGN ACCURACY

- Single tubing with connectic with screw
- "Watchmaking" mechanism

#### SHOCK RESISTANT

- Non deformable metal casing
- Shock proof foam manometer protector

#### COMFORT

- Ambidextrous thanks to the central decompression valve
- Large pump for quick inflation
- SPENGLER Velcro fastening

#### CUFF

- Coton or or nylon cuff with velcro fastening
- Sizes: XS (10-15 cm), XS (15-22 cm), S (20-28 cm), M (26-34 cm), L (32-42 cm) ou XL (42-52 cm)
- Available in blue, grey and Cartoon
- Available in Clinic® to prevent nosocomial infection
- Comes with a transport pouch

Adult (M), cotton, velcro, navy: ref 520627

Adult (M), cotton, velcro, grey: ref 520621

Adult (L), nylon, velcro, grey: ref 520620

Child/Small adult (S), cotton, velcro, motifs: ref 520622

Multicuff (L,M,S), velco: ref 520630

### Features

|                   |               |                |                               |
|-------------------|---------------|----------------|-------------------------------|
| CE Marking        | CE0459        | Medical device | Class Im (measuring function) |
| Bar code          | 3700446008864 | HS Code        | 90189010                      |
| Warranty (months) | 36            | Country        | France                        |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.691 |

#### SPECIFIC STORAGE CONDITIONS

|                     |     |                 |    |
|---------------------|-----|-----------------|----|
| Minimum temperature | -20 | Max temperature | 70 |
|---------------------|-----|-----------------|----|

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)