

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. D1.212.212



 12/08/2025 21:49:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

LUXAMED®

Ref. D1.212.212

Penlight for eye and throat diagnostic., long life technology.

Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

- White: D1.212.912
- Blue: D1.212.212
- Grey: D1.212.312
- Yellow: D1.212.512
- Kids : D1.312.912
- Black: D D1.212.112
- Red : D1.212.412
- Green : D1.212.612

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961758	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.04		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			

Colors: Blue
08/2025 21:49:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr