

ECG PATIENT CABLE 4.5M (FOR TANGO M2 ONLY)

SUNTECH

Ref. DSU230717



06/08/2025 22:20:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ECG PATIENT CABLE 4.5M (FOR TANGO M2 ONLY)

SUNTECH

Ref. DSU230717

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012304	HS Code	9018110000
Warranty (months)	3		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


06/08/2025 22:20:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr