

DISPOSABLE SPECULUM DIA. 4 MM, BOX OF 250, BLACK, ADULT, SPENGLER



Ref. L20748



06/08/2025 22:18:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DISPOSABLE SPECULUM DIA. 4 MM, BOX OF 250, BLACK, ADULT, SPENGLER



Ref. L20748

Description

Box of 250.
For adult, single use,
Black,
Ø 4 mm.
Universal, compatible with all brands
Box with dispenser spout for easy access to speculums

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446031275	HS Code	3926909790
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	11
Packed unit / Length (cm)	11	Packed unit / Weight	0.135
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/08/2025 22:18:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr