

BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-250



05/08/2025 23:05:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

BABINSKY HAMMER, BLACK COLOR WITH PLASTIC HANDLE S.S



Ref. M05-320-250

=====

Description

- Handle with tip on the screw
- Unscrewable head 2 positions
- Length 25 cm
- Available in blue (M05-320-260), black (M05-320-250) and red (M05-320-280).

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446010652	HS Code	90189090
Warranty (months)	60	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	25
Packed unit / Weight	0.089	Bare Unit / Weight (kg)	0.089
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

📅 05/08/2025 23:05:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
🌐 www.gsh-med.fr

📞 +33(0)4 42 90 31 31
✉ info@gsh-med.fr