

Ref. M21000



06/08/2025 22:32:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

Ref. M21000

Description

Pressure cuffs are medical devices designed for use in operating theatres, intensive care units and emergency transport vehicles. They are used to manually accelerate the infusion of blood and solutions when gravity is not sufficient to achieve the desired flow rates.

The GREY FUSE pressure cuff is compatible with all IV and blood bags, and comes in 3 sizes: 500cc (ref. M20500) , 1000cc (ref. M21000) and 3000cc (ref. M25000).

Optimised monitoring and reliability

Reusable and easy to clean/disinfect

Visualisation of the IV bag thanks to a transparent fabric made up of a screen and a polyamide mesh.

Pressure range from 0 to 300mmHg in 10mmHg increments

Pressure gauge accuracy ± 3 mmHg

Red zone on dial at 300mmHg to prevent overpressure

Pressure relief controlled by a simple tap.

All our cuff are made in our Castelsulnier factory near Agen.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468211501	HS Code	9018908400
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	470		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40
OTHER FEATURES			

2-28 - For information purposes only, can be subject to change