


## NC3 MONITOR NIBP+SPO2 COMEN+TEMP

Ref. BCMONI02

# COMEN



CE  
1639



DISPOSITIF MEDICAL  
CLASSE  
**IIb**  
MEDICAL DEVICE

GARANTIE  
3 ANS  
**3**  
GUARANTEE



12/08/2025 22:00:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

# NC3 MONITOR NIBP+SPO2 COMEN+TEMP

# COMEN

Ref. BCMONI02

**The compact, fast and accurate diagnostic station, perfect for outpatient bedside with temperature option**

## Description


The NC3 monitor measures 3 vital parameters: non-invasive blood pressure (NIBP), oxygen saturation (SpO2) and temperature (TEMP) using an Infra Red ear thermometer. Particularly fast and easy to use (starts in 1 to 2 seconds, PNI measurement: 20 to 30 seconds) Its weight and small size with optional trolley with wheels or a wall mount make it a very suitable instrument for outpatient at the bedside, in hospitals, clinics or nursing homes, as well as in doctor's private practices.


- 6" LED backlit screen
- PNI and SpO2 accessories for all types of patients (adult, pediatric, neonatal) (delivered with adult PNI cuff and adult SpO2 sensor)
- IR thermometer
- Measurement of spot check (manual) or programmed (automatic) NIBP
- Integrated rechargeable lithium battery for a continuous autonomy of 12 hours
- Memory of 50 groups / patients data
- Weight: 1.25kg
- Options: Trolley, wall mount

## Features

CE Marking	CE1639	Medical device	Class IIb
Bar code	3661040219066	HS Code	9018191000
Warranty (months)	36	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	12.5
Packed unit / Height (cm)	26	Bare Unit / Height (cm)	22
Packed unit / Length (cm)	31	Bare Unit / Length (cm)	13
Packed unit / Weight	2.95	Bare Unit / Weight (kg)	1.25
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	55
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 12/08/2025 22:00:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)