


BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

Ref. MDOUES01

De Boissy
— DEPUIS 1982



06/08/2025 22:28:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BAG OUESSANT, DE BOISSY, WITH BLACK LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDOUES01

Eco-responsible, customizable and "super organized" for a trendy look with all the practical and modern options

Description

All fabrics, exterior and interior, certified Global Recycled Standard (GRS), have been made from recycled material and remain washable, resistant and water-repellent. Its deep blue, inspired by the shores of Ushant, is enhanced by black grained leather-effect finishes for an elegant, professional look. If De Boissy wanted to pay tribute to the caregivers and be inspired for his bags in the shape of the traditional leather bag of doctors, the look of Ouessant bags contrasts well and truly with the usual classicism. De Boissy played the card of elegance, trend and modernity. And the fact remains that they have been designed to be "ultra practical", to be the best ally of health professionals during their medical visits at home or in community. - The large volume interior main space includes 6 pockets, 2 removable velcro dividers, 1 padded compartment for laptop. With a very large bidirectional zip opening, it gives very easy access to diagnostic instruments and care equipment - This space also includes 2 detachable compartments (with velcro): 1 transparent to quickly identify the contents, and 1 isothermal, supplied with its reusable cold gel pocket, to slide up to 20 bulbs. There is also a carabiner to, for example, gather the keys of patients to visit every day and even an elastic buckle to hold a bottle or a gourd upright. - The ventral pocket (40x20 cm) opens completely with on one side 2 zipped storage nets and 2 flat pockets, and on the other side several elastic loops allow to store small care supplies: pliers, scissors, pen lamp, thermometer, bandages, pens, etc ... - Three external zipped pockets, one for A4 documents and another on the side with elastic to slide the bottle of hydroalcoholic gel which remains easily accessible at any time - 1 pouch, grained leather appearance, is detachable. What allows you to customize your bag (fancy printed pouch sold separately) - An external USB port with integrated charging cable allows you to connect a backup battery and never fail your phone, or forget your charger at a patient - Padded and solid handles, a padded and adjustable shoulder strap - A name tag to identify the bag - A non-slip bottom to reduce friction and wear - Dimensions w x h x d: 40 x 35 x 24 cm - Weight: 1.6 kg - Maximum recommended content: 6 kg / 24 liters - Sold empty The De Boissy bag is put back in its dust bag to protect it when stored. Grey with patterns, this accessory is designed from recycled materials with a cord tightening.

📅 06/08/2025 22:28:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDOUES01

Features

CE Marking	NONE	Medical device	No
Bar code	3700446042417	HS Code	42022290
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	46	Packed unit / Height (cm)	44
Packed unit / Length (cm)	60	Packed unit / Weight	1.65
Bare Unit / Weight (kg)	1.65		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

06/08/2025 22:28:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr