

## ISOTHERMAL BAG ELITE BAG MINICOOL'S, BLUE

**ELITE BAGS**

Ref. MECOOL04



05/08/2025 22:51:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ISOTHERMAL BAG ELITE BAG MINICOOL'S, BLUE



Ref. MECOOL04

### Description

Insulated bag that combines the portability and practicality of ROW'S, with the capacity of COOL'S

The bag is a perfect size for laboratories transporting medium quantities of samples at home, as well as for samples from hospitalized patients, and even for vaccination campaigns.

- 1 isothermal main compartment with 4 terry supports to accommodate 44 tubes of different diameters, with 4 hermetic pockets to prevent liquid leakage.
- 1

insulated compartment for vials, with removable separators that prevent spilling•

The insulated compartments have pockets to place cold accumulators• 1

compartment with transparent front side for laboratory equipment•

1 external digital thermometer to monitor the temperature inside•

Outside, a transparent pocket for documents•

Stackable bag

- Material: 420D Polyester

- Outdoor measures: 45 x 24 x 27 cm / 29.2L•

Weight: 0.48 kg

- Recommended max. capacity: 12 kg

### Features

CE Marking	NONE	Medical device	No
Bar code	3661040167459	Warranty (months)	24
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

05/08/2025 22:51:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)