

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_L



07/08/2025 22:24:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

PATIENT STOOL WITH NYLON BASE 5 DOUBLE CASTORS BLUE



Ref. OS007_L

Plastic base stool (seat 3.5 cm)

Description

- Seat ø: 43 cm, padded in very high density foam (30 kg / m3)
- Fireproof seat (fire resistance classification M1)
- Axis, tube ø: 5.5 cm in stainless steel
- Black plastic base with 5 branches ø: 58 cm, very stable with double casters
- Height adjustable by simple pressure and gas piston
- Height variation from 53 to 70 cm
- Maximum load: 150kg

Matching stool to the OSA003 sofa (several colors available for both models)

- Available in Black (OS007_N) and Lavender (ref. OS007_L)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040196138	HS Code	9402900000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Packed unit / Height (cm)	42
Packed unit / Length (cm)	41	Packed unit / Weight	6.2
Bare Unit / Weight (kg)	5.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Lavender	Packed unit / Quantity	1

07/08/2025 22:24:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr