

LUXASCOPE AURIS SET LED 3.7 V, WHITE (OTOSCOPE + OPHTHALMOSCOPE)

LUXAMED®

Ref. AF.426.914

**CE**

 14/08/2025 23:31:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUXASCOPE AURIS SET LED 3.7 V, WHITE (OTOSCOPE + OPHTHALMOSCOPE)

LUXAMED®

Ref. AF.426.914

Complete rechargeable set for precision diagnosis: Auris LED 3.7V black otoscope + LuxaScope LED 3.7V black ophthalmoscope

Description

Otoscope :

Lighting: LED ring with 6 LEDs (concentric),

no fibre optics, the field of vision is not limited.

- Maintenance-free, as no bulb change is required

- Possibility of varying brightness intensity

- Automatic switch-off after 3 minutes

- LED life: approx. 100,000 hours at constant current

- Illuminance: approx. 10,000 Lux

- Colour temperature: approx. 4,200 K

- Optimised head design

- Swivelling magnifier, 3x magnification

- Injection port for pneumatic test

- Handle made from two materials: aluminium and fibreglass

reinforced, compatible with dermatoscope head and ophthalmoscope

LuxaScope LED

- ON/OFF push-button switch

- Power supply: 2x AAA alkaline batteries, included

- Includes 20 grey disposable specula (10 x 2.5mm and 10 x 4.0mm)

Ophtalmoscope :

Correction values: + 20 bis - 20 D (± 1 2 3 4 6 8 10 15 20 D)

- 5 different filters available: large spot, small spot, green filter, half-moon and reticle

- Illumination: LED

Brightness can be dimmed to avoid glare

14/08/2025 23:31:09 - For information purposes only, can be subject to change

- Automatic switch-off after 3 minutes

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

reinforced, compatible with otoscope and dermatoscope heads

LUXASCOPE AURIS SET LED 3.7 V, WHITE (OTOSCOPE + OPHTHALMOSCOPE)

LUXAMED®**Ref. AF.426.914**

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964513	HS Code	9018908400
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

 14/08/2025 23:31:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr