

## CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM

Ref. DUCRYO03



14/08/2025 23:05:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CRYO PROFESSIONAL 170 ML - 30 APPLICATIONS 5 MM & 2 MM



Ref. DUCRY003

Effectively treats benign skin lesions

### Description

- Cryotherapeutic agent, allowing to deliver a temperature of -55 ° C in less than 15 seconds
- Application time between 15 and 40 seconds
- Complete set including:
  - A 170 ml aerosol can
  - 2mm and / or 5 mm applicator sticks, usable depending on the different lesions
  - An instruction manual
- The sale of this medical equipment is reserved for health institutions, doctors and doctors' offices

### Features

CE Marking	CE0344	Medical device	Class IIa
Bar code	8717484005835	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	14
Packed unit / Length (cm)	14	Packed unit / Weight	0.3
SPECIFIC STORAGE CONDITIONS			
Max temperature	25		
OTHER FEATURES			

14/08/2025 23:05:27 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)