

PHIAL AMPOULARIUM ELITE BAGS EB09.002



Ref. MEAMPO01



14/08/2025 23:33:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

PHIAL AMPOULARIUM ELITE BAGS

EB09.002



Ref. MEAMPO01

Insulated amp maker, 43 bulbs

Description

- Padded central separator
-

Fasteners- Color: black

- Weight: 160 g
-

Outdoor measures: 18 x 13 x 4 cm / 0.93 L

- Maximum recommended capacity: 1 kg

- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540546	HS Code	4202199090
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	13
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	18	Bare Unit / Length (cm)	18
Packed unit / Weight	0.276	Bare Unit / Weight (kg)	0.16
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 14/08/2025 23:33:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr