

COOL BAG FOR DIABETIQUE PM DIA, GREY

 **ELITE BAGS**

Ref. MEDIAB03_B



 14/08/2025 23:38:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

COOL BAG FOR DIABETIQUE PM DIA, GREY



Ref. **MEDIAB03_B**

Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner

Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening
- Material: Polyester graybitone
- Weight: 0.12 kg.
- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03_N), blue jeans (MEDIAB03_B) or pink (MEDIAB03_RC)
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542458	HS Code	4202.12.99.90
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	17
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.012
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Bleu jean	Packed unit / Quantity	1

14/08/2025 23:38:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr