

LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

EMERGENCY'S™

Ref. MEEMER03_R



14/08/2025 23:24:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

EMERGENCY'S™

Ref. MEEMER03_R

Intended to transport the most fragile + relief materials

Description

- Four compartments for interior supports and high-capacity compartment.

- Transparent pouch in the inner flap, designed to see and organize the material easily, thanks to the rubber-elastic that guarantees a correct fixation- Five outer pockets,

two side pockets, one on the back and one on the top and a front pocket that can be divided internally thanks to rubber pockets and bands

- Base of the waterproof bag

- High resistance zips

- Reflective

strips- Reinforced handles and shoulder straps

- External card holder

- Dimensions: 55 x 31.5 x 33 cm, capacity 38 L

- Internal dimensions: 40 x 20 x 24 cm

- Material 600D Polyester

- Weight 1.76 kg

- Sold empty

- Available in : Red (MEEMER03_R), Blue (MEEMER03_B) and Orange (MEEMER03_O)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540799	HS Code	420299000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	28	Bare Unit / Width (cm)	31.5
Packed unit / Height (cm)	20	Bare Unit / Height (cm)	33
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	55
Packed unit / Weight	1.8	Bare Unit / Weight (kg)	1.76
SPECIFIC STORAGE CONDITIONS			

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr