


WALL UNIT WITH OTOSCOPE-MODULE, OPHTHALMOSCOPE-MODULE AND EAR FUNNEL

LUXAMED®

Ref. W1.336.114



 14/08/2025 23:32:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WALL UNIT WITH OTOSCOPE-MODULE, OPHTHALMOSCOPE-MODULE AND EAR FUNNEL

LUXAMED®

Ref. **W1.336.114**

Wall Unit with Otoscope-module, Ophthalmoscope-module and ear funnel

Description

The wall station is equipped to provide the required diagnostic aid, which is useful in daily life in a hospital, clinic, or private practice.

- Resealable dispenser of speculums with a lid and can hold up to 30 speculums per dispenser.
- Power supply on mains 220V / 110V : 1 power cable + 1 EU/GB/US adapter included
- Automatic shut-off after 3 min of use
- Possibility to vary the brightness with a simple long press

Available in :

- Wall station with otoscope and speculum dispenser*
- Wall station with a headless sleeve

*Speculums not included with the dispenser

Features

CE Marking	NONE	Medical device	No
Bar code	4260381965008	HS Code	940549400
Country	Germany		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

14/08/2025 23:32:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr