


WALL UNIT MODULE WITH 1 HANDLE ONLY WITHOUT INSTRUMENT HEAD

LUXAMED®

Ref. W9.316.114



 14/08/2025 23:55:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WALL UNIT MODULE WITH 1 HANDLE ONLY WITHOUT INSTRUMENT HEAD

LUXAMED®

Ref. **W9.316.114**

Wall Unit module with 1 handle only without instrument head

Description

The wall station is equipped to provide the required diagnostic aid, which is useful in daily life in a hospital, clinic, or private practice.

- Power supply on mains 220V / 110V : 1 power cable + 1 EU/GB/US adapter included
- Automatic shut-off after 3 min of use
- Possibility to vary the brightness with a simple long press

Available in :

- Wall station with otoscope, ophthalmoscope and speculum dispenser*
- Wall station with otoscope and speculum dispenser*

*Speculums not included with the dispenser

The device is sold with the handle only, and must be completed with an otoscope, dermatoscope or ophthalmoscope head.

Features

CE Marking	NONE	Medical device	No
Bar code	4260381965022	HS Code	940549400
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

📅 14/08/2025 23:55:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr