


OSCAR 2 ABP MONITOR - STANDARD

Ref. DSU230220

SUNTECH



CE
0413



OSCAR 2 ABP MONITOR - STANDARD

Ref. [DSU230220](#)

SUNTECH

The reference for ambulatory blood pressure monitoring

Description

- Clinically tested by the British Hypertension Society and the European Society of Hypertension, the Oscar 2 M250 offers very high precision measurements- SphygmoCor technology (optional) that provides important information on the clinical impact of arterial stiffness and pulse wave reflection, allowing an informed examination of cardiovascular risk
- Extraordinary comfort thanks to the ORBIT cuff
- Silent
- Lightweight and compact
- A mobile application so that the patient can inform his activity via Bluetooth (optional)
- Up to 48 hours of recording
- Analysis and editing software Accuwin Pro 4.0

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	0110840935105978	HS Code	90181990
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Packed unit / Height (cm)	19
Packed unit / Length (cm)	30	Packed unit / Weight	1.815
Bare Unit / Weight (kg)	0.284		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--