

CABLE PATIENT ECG TANGO

Ref. DSU230716

SUNTECH

Cable patient ECG Tango

Description

ECG cable for Tango® automatic monitor

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012298	HS Code	9018110000
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.217

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---