

ECG PATIENT CABLE 4.5M (FOR TANGO M2 ONLY)

SUNTECH

Ref. [DSU230717](#)

ECG patient cable 4.5m (for Tango M2 only)

Description

ECG patient cable 4.5m (for Tango M2 only)

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012304	HS Code	9018110000
Warranty (months)	3		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--