

Ref. E1.416.012



CE



Description

07/10/2025 22:55:09 – For information purposes only, can be subject to change

Ref. E1.416.012

Features

CE Marking	CE	Medical device	Class I
Bar code	4260381960331	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES