

EMERGENCY'S RESCUE BACKPACK - RED WATERPROOF

Ref. MEEMER05_WA

EMERGENCY'S™



07/10/2025 23:54:41 – For information purposes only, can be subject to change

EMERGENCY'S RESCUE BACKPACK - RED WATERPROOF

EMERGENCY'S™

Ref. MEEMER05_WA

Its foldable flap opening provides access to the tools necessary for emergency actions

Description

- Main ergonomic handle with Velcro- Large pockets on both sides
- Padded shoulder straps and backs, adjustable backpack
- Large front pocket, ready to accommodate a defibrillator
- Spacious interior for convenient and easy access to the material
- Internal compartments with mesh and transparent windows for the organization of the material
- External dimensions: 37 x 45 x 21 cm / 30L
- Material: 600D Polyester and PVC
- Weight: 1.22 kg
- Sold empty

EMERGENCY'S RESCUE BACKPACK - RED WATERPROOF

EMERGENCY'S™

Ref. [MEEMER05_WA](#)

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542748	HS Code	4202929890
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	32	Bare Unit / Width (cm)	37
Packed unit / Height (cm)	33	Bare Unit / Height (cm)	21
Packed unit / Length (cm)	55	Bare Unit / Length (cm)	45
Packed unit / Weight	1.76	Bare Unit / Weight (kg)	1.22

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---