

CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS

Ref. 260012



CE
0123



08/10/2025 02:20:26 – For information purposes only, can be subject to change

CARDIOMATE 12 - ELECTROCARDIOGRAPH 12 CHANNELS



Ref. 260012

CARDIOMATE® ECGs are highly reliable instruments.

Description

ECG Cardiometre 12® are robust and reliable ECG. Featuring a 10" color screen tilt they allow viewing 12 lead simultaneously. Simple to use, recording and printing ECG are triggered with only one button. The Cardiometre 12® are Cardiometre equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer. The Cardiometre 12® are equipped with automatic interpretation and pacemaker detection. The internal battery life of 3 hours, their lightness and carrying handle make Cardiometre 12® the perfect companion in the office or traveling.- Carrying handheld.- Tilting screen :10-inch color LCD display (resolution: 800x480 px).- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM.- Autonomy :The operating life is 1080 ECG.- Automatic interpretation and intervals calculation.- Pacemaker detection.- Network connection, USB & printer.- Printing and rapid analysis.- Acquisitions and display 12 simultaneous tracks.- Alphabetic keypad.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446025496	HS Code	9018110000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	40
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	9.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	42
Packed unit / Weight	8.3	Bare Unit / Weight (kg)	4.2

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	50
Minimum humidity	25	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---