

WIFI ANTENNA USB FOR ECG CARDIOMATE SPENGLER



Ref. 260013

Wifi Antenna USB for ECG Cardiomate SPENGLER

Description

USB WIFI antenna for ECG Cardiomate 3, 6 and 12 tracks.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446042080	HS Code	9018110000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.005		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--