

Ref. 260020

08/10/2025 01:22:07 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

**PATIENT CABLE FOR ECG CARDIOMATE *
SPENGLER**Ref. 260020

Patient Cable for ECG CARDIOMATE * SPENGLER

Description

Patient Cable for ECG CARDIOMATE * SPENGLER

Features

CE Marking	NONE	Medical device	No
Bar code	3700446025502	HS Code	8544200090
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.4

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---