

## LIMB ELECTRODES CLAMP FOR ECG CARDIOMATE X 4

Ref. 260021



LIMB ELECTRODES CLAMP FOR ECG CARDIOMATE x 4

### Description

Kit of 4 limb clips for ECG Cardiomate

08/10/2025 01:39:38 – For information purposes only, can be subject to change

## LIMB ELECTRODES CLAMP FOR ECG CARDIOMATE X 4

Ref. 260021



### Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 3700446025519 | HS Code        | 9018110000 |
| Warranty (months) | 60            | Country        | China      |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |     |
|---------------------------|----|---------------------------|-----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc  |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.2 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|