

BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON

Ref. 260023



BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON

Description

Kit of 10 adapters for snap ECG electrodes CARDIOMATE

08/10/2025 02:22:22 – For information purposes only, can be subject to change

BAG OF 10 UNITS ECG ADAPTER PRESSURE BUTTON

Ref. 260023



Features

CE Marking	NONE	Medical device	No
Bar code	3700446030704	HS Code	9033009000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.06

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---