

3-CHANNEL ECG PAPER, ROLLED

Ref. 260024



3-CHANNEL ECG PAPER, ROLLEDRef. 260024

=====

Description

PAPER ROLL dimensions 80MM x 20M FOR CARDIOMATE ECG 3

Features

CE Marking	NONE	Medical device	No
Bar code	3700446026677	HS Code	4823400000
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	8	Packed unit / Weight	0.102

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---