

6-CHANNEL ECG PAPER, ROLLED

Ref. 260025



08/10/2025 00:37:30 – For information purposes only, can be subject to change

6-CHANNEL ECG PAPER, ROLLED

Ref. 260025



6-channel ECG Paper, Rolled

Description

Paper roll dimensions 110mm x 20m for ECG Cardiomate 6

Features

CE Marking	NONE	Medical device	No
Bar code	3700446027377	HS Code	4823400000
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	12	Packed unit / Weight	0.134

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---