


CARDIOMATE EVI 7" TABLET ECG SPENGLER

Ref. 260050



CE
0123



07/10/2025 23:44:44 – For information purposes only, can be subject to change

CARDIOMATE EVI 7" TABLET ECG SPENGLER

Ref. 260050



Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 7" tablet, a fully autonomous technological product.

Description

Very easy to use:

- 7" color touch screen
- Portable: compact size and light weight, easy to carry: 0.8kg and 197 x 112.4 x 26.1 mm
- Large battery capacity, rechargeable on mains
- Allows to print ECGs directly on a printer via wifi network

High-performance technology:

- High accuracy: sampling rate 32,000 Hz
- Detection and marking of the pacemaker
- Automatic interpretation (Glasgow) and calculation of intervals
- Input of patient information using a keyboard or barcode reader
- real time visualization of the 12 tracks at the same time as automatic interpretation and interval calculations. (possibility of a 18 leads tablet optional)
- Data export via wifi or USB port under .pdf or .xml

CARDIOMATE EVI 7" TABLET ECG SPENGLER

Ref. 260050



Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6971049341789	HS Code	9018110000
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Packed unit / Height (cm)	28
Packed unit / Length (cm)	33	Packed unit / Weight	2.8
Bare Unit / Weight (kg)	0.8		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES