



SPENGLER STETHOSCOPE PULSE II SINGLE, PAEDIATRIC

Ref. 507704



07/10/2025 22:50:06 – For information purposes only, can be subject to change

**SPENGLER STETHOSCOPE PULSE II
SINGLE, PAEDIATRIC****Ref. 507704**

Stethoscope, binaural, complete: A mechanical listening device designed for listening to sounds from the heart and lungs from paediatric patients.

Description

HIGH ACOUSTIC QUALITY

- Aluminium single head chestpiece
- High performance single tubing
- Aluminium binaural
- High performance diaphragm

COMFORT

- Non-chill ring
- Silicone earpieces
- New binaural design better adapted to human anatomy

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029906	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	3
Packed unit / Length (cm)	24	Packed unit / Weight	0.215
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		