


## LIAN METAL / GREY NYLON CUFF / OBESE (XL)

Ref. 521317



CE  
0459



08/10/2025 00:16:03 – For information purposes only, can be subject to change

**LIAN METAL / GREY NYLON CUFF / OBESE  
(XL)****Ref. 521317**

=====

## Description

Large dial, ambidextrous model (simply rotate to change hands). Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

Available in :

Adult (M), cotton, velcro, red (521317)  
Adult (M), cotton, velcro, grey (521311)  
Adult large (L), cotton, velcro, grey (521310)  
Child/Small Adult (S), nylon, velcro, grey (521316)  
Multibrassard (S, M, L) velcro (521330)

## Features

|                             |               |                           |                               |
|-----------------------------|---------------|---------------------------|-------------------------------|
| CE Marking                  | CE0459        | Medical device            | Class Im (measuring function) |
| Bar code                    | 3700446008857 | HS Code                   | 90189010                      |
| Warranty (months)           | 36            | Country                   | France                        |
| DIMENSIONS AND WEIGHT       |               |                           |                               |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc                            |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.75                          |
| SPECIFIC STORAGE CONDITIONS |               |                           |                               |
| Minimum temperature         | -20           | Max temperature           | 70                            |
| Minimum humidity            | 85            | Maximum humidity          | 85                            |
| OTHER FEATURES              |               |                           |                               |
| Packed unit / Quantity      | 1             |                           |                               |