

MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



CE
0123



07/10/2025 23:07:18 – For information purposes only, can be subject to change

MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.

Description

HIGH MEASUREMENT

ACCURACY• Double tubing for linear pressure rise

COMFORT• Ambidextrous•

Remotely readable measurement•

2.5-meter spiral tubing• Large volume ergonomic pear for quick inflation•

Storage basket behind the dial• Comes with cuff wheeled base

- Comes with armband size M, nylon

OPTION

- 3 sizes of cuffs (S, M, L) available

MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446028152	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.16
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		