


Ref. C1.616.114

LUXAMED®



Technologie CCT
2 couleurs de lumière



08/10/2025 01:44:44 – For information purposes only, can be subject to change

=====

Ref. C1.616.114

LUXAMED®

 08/10/2025 01:44:44 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

Ref. C1.616.114

Innovative dermatoscope in pocket size

Description

Innovative illumination: LED-Ring with 2x5 LEDs in concentric structure for a highly detailed image, homogenous light, switchable with one button press : warm white light or cold white light

- Removable contact plate for easy cleaning
 - Maintenance-free, because no lamp replacement necessary
 - Stepless zoom area with 10 x magnification for a distortion-free image
 - Optimized design of head
 - Automatic switchoff after 3 minutes
 - Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
 - Intensity of illumination: 10.000 Lux
 - Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
 - Two-component-handle (aluminium/fibre-glass reinforced),
 - Also combinable with LuxaScope Auris CCT LED 2.5 V
- otoscope-head
- Electronic ON/OFF push-button-switch
 - Power supply: 2 x AAA Alkaline batteries, inclusive

Available in :

- Black : C1.616.114
- White : C1.616.914

=====

Ref. C1.616.114

LUXAMED®

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963615	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.14		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Black		

08/10/2025 01:44:44 – For information purposes only, can be subject to change