

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. [DS99406](#)



Description

Adult cuff (26 to 35 cm)

Washable nylon cuff without pocket with Velcro fastening

08/10/2025 02:46:36 – For information purposes only, can be subject to change

ADULT CUFF COMPLETE WITH QUICK CONNECT FOR EASY2/EXPERT(26 TO 35 CM)



Ref. [DS99406](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040075211	HS Code	90189010
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.058	Bare Unit / Weight (kg)	0.058

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---